

# Forensic analyzer systems



*multi-purpose laboratory analyzer*



*portable analyzer with remote fiber-optic sensing*



*vehicle mounted system*

A unique range of Raman spectroscopy systems designed for use by law enforcement agencies

- multipurpose Raman analyzer for laboratory use
- portable Raman analyzer for on-site/scene of crime use
- advanced UV Raman systems for specialist applications

## **Applications**

- narcotic identification
- explosive identification
- fiber identification
- paints/pigments analysis
- gunshot residue analysis
- gemstone identification

## **Advantages**

- non-destructive analysis
- analysis through transparent packages, including evidence bags
- remote optical fiber analysis (over 100 m)
- fast material identification with extensive software libraries
- no sample preparation needed
- samples as small as 1  $\mu\text{m}$  can be studied

## **Usage**

- forensic applications worldwide

systems for law enforcement agency use

**RENISHAW** 

# Forensic analysis

## The advantages of the technique...

- state-of-the-art Raman technology
- unique characterization of materials
- rapid analysis
- non-destructive - evidence is preserved
- no sample preparation required
- material inside most transparent packaging can be analyzed



## ...and of the analyzers

- easy-to-use automated systems
- rugged
- fast material identification using extensive software libraries
- traceability and data transfer between different instruments
- in-built automated calibration routines
- high spatial resolution (to 1  $\mu\text{m}$  for laboratory analyzer)
- spectral libraries for material identification

## in the laboratory...



- easy sample loading - as simple as using an ordinary optical microscope
- video viewer for rapid positioning of the sample, with image capture facility
- automated sample stage for 2-D mapping of surfaces with focus-tracking for rough surfaces

## ...and in the field



### Customization options:

- mobile operation (trucks etc.)
- extra rugged for in-field use
- stabilized mounts for fiber optic probes
- fiber probes with remote viewer for probe positioning
- battery operation developments

## Raman spectroscopy

### Complementary techniques: Raman and IR

Raman scattering and infrared absorption are complementary techniques: both use light to probe the molecular structure of materials, and can identify the composition of unknown substances.

In infrared absorption, light is absorbed when its frequency matches that of one of the bonds in the molecule. This absorption occurs at relatively long wavelengths, in the infrared region of the electromagnetic spectrum.

In the Raman case, the molecule scatters incident laser light, reducing its frequency as bonds are set vibrating. The reduction in frequency of the light corresponds to the vibration frequency of the bond.

The majority of lasers used for Raman spectroscopy have visible or near-visible emission frequencies, resulting in scattered light frequencies that are also visible or near-visible. Instruments using ultraviolet excitation are also available for specialist applications.



### Differences between Raman and IR

Raman has advantages over infrared absorption

- It is non-destructive
  - samples require no preparation
  - technique is non-contacting
  - evidence is preserved
- High temperature studies can be performed easily
- Lower Raman shift coverage; 100  $\text{cm}^{-1}$  possible with standard optics. This allows a wider range of materials to be studied
- Microscopy:
  - smaller particles can be studied than IR (1  $\mu\text{m}$  to 2  $\mu\text{m}$  across, an order of magnitude better than IR)
  - confocal measurements are easy (high depth and lateral resolutions) allowing embedded particles to be analyzed
- Sampling is easier:
  - water is the ideal solvent, glass the ideal container/window material
  - samples within transparent plastic containers can be studied
  - a wide variety of sample shapes and sizes can be accommodated

### Raman and IR





Diagram illustrating infrared absorption and Raman scattering.

### Raman and IR spectra



Raman and infrared spectra of PETN explosive, showing the correspondence between the two.

Data courtesy of Prof. D.N. Batchelder, University of Leeds.

# Forensic analysis

## Narcotics/drugs

The Renishaw forensic analyzers are ideal for detecting and identifying all the common narcotics, including:


- cocaine, and cutting agents (sugar, chalk, etc.)
- amphetamines; characterization of manufacturer source

Narcotics can be identified even after they have been “cut” with other substances.



### Detection


combination of Raman images taken at 1601  $\text{cm}^{-1}$  and 1713  $\text{cm}^{-1}$



The Renishaw laboratory analyzer has a unique filter system that allows fast direct 2D imaging. This facility can be used to locate rapidly the suspect material.

Images of a mixture of cocaine hydrochloride and sugar are shown above. The standard microscope image (left) shows a mass of indistinguishable material, whereas the image taken with the laboratory analyzer (right) reveals particles of cocaine as bright regions.


### Identification



Spectra of amphetamines.

Data courtesy of Profs. Zhao JT and Zhang PX,  
Yunnan Polytechnic University, Kunming, PR China.

### Identification




Spectra of morphine and heroin.

Data courtesy of Profs. Zhao JT and Zhang PX,  
Yunnan Polytechnic University, Kunming, PR China.

The Renishaw Raman forensic analyzers can analyze all of the common explosives and associated materials (TNT, Semtex, ammonium nitrate/sugar, etc.)

Traces of explosive can be detected non-destructively and identified on a wide range of materials, including fabrics, metals and paintwork, without the need for sample preparation.


### Identification




Raman spectra of PETN and RDX explosives (each from approximately 10 pg of material). The spectra are unique to the material, allowing explosives to be distinguished from each other, and from other materials.

### Detection

white light




Raman image using strong RDX band



Images of a fingerprint on aluminum taken with a standard microscope (left), and with the laboratory analyzer (right). A particle of RDX explosive, a few micrometers across, is clearly revealed on the Raman image from the analyzer.

### Remote detection



Raman spectra of PETN explosive.

Top: bulk material spectrum taken with the laboratory analyzer.  
Bottom: traces of explosive in a fingerprint examined using a remote fiber optic probe (finger grease features are also present in the spectrum)


# Forensic analysis

## In situ studies through packaging

Renishaw's Raman systems can acquire spectra from samples contained within a wide range of transparent/translucent packaging materials, such as glass vials and plastic bags.

This allows materials (such as narcotics contained within evidence bags) to be identified without disturbance.

### Identification



Spectrum of cocaine taken through a plastic (polyethylene) bag. The Raman bands of the cocaine appear on a luminescent background from the bag (633 nm excitation).




## Paints and pigments


Spectra taken from paint particles differ widely, depending on the paint type and manufacturer. A typical application of this capability is the linking of tiny paint particles to a vehicle involved in a road accident.

The high spatial resolution of the laboratory analyzer allows different layers of the paint fragment (corrosion inhibitor, undercoat, etc.) to be studied, giving even greater discrimination.

### Identification



Typical paint spectra, in this case from three different green paints. They can easily be distinguished from each other.




Data courtesy of Niu TL, Police Ministry, PR China.

## Forgery of printed material (banknotes etc.)

The advent of inexpensive printing and publishing technology has enabled criminals to forge convincing documents with great ease.


The Renishaw forensic analyzers can be used to distinguish between pigments, paints, and papers that appear identical. This can all be performed rapidly, without damaging the documents.

### Identification




Spectra of inks on a \$10 (US) banknote.

### Identification



Spectra of inks on a 100 Yuan (Chinese) banknote.

### Identification




Spectra of inks on a £10 (Sterling) banknote.

# Forensic analysis

## Handwriting forgery

Many inks appear identical when viewed with the naked eye, or with a conventional microscope, making it difficult to detect forged signatures and other such document adulteration.


However, the molecular compositions of inks vary markedly, allowing the Renishaw forensic analyzers to distinguish between them. The high spatial resolution of the Renishaw analyzers enables inks to be identified on paper and other materials.



*Data courtesy of Niu TL, Police Ministry, PR China.*

## Fiber identification

The high spatial resolution of the laboratory analyzer makes it ideally suited to identifying single fibers that have been transferred during crimes. The resulting spectra enable not only the fiber material but also the fiber pigments to be identified.



## Synthetic minerals and gems


The Renishaw forensic analyzers are ideal tools for identifying gemstones.

Renishaw's **inorganic database** of spectra enables the composition of suspect gemstones to be determined rapidly.

The high resolution of the laboratory analyzer's microscope also allows users to detect and identify microscopic flaws in gemstones, and the resins and fillers used to hide them. The same technique can be used to study inclusions in gems to determine their source.



### Detection of gem flaws




Spectrum taken from a faint feature in a diamond (inset image) reveals that a resin has been used to fill a crack in an attempt to hide it.


*Data courtesy of J Yuan International Inc., Taiwan.*

## Gunshot residues

Gunshot residues consist of materials such as carbonates, nitrates, and oxides, (originating from the explosive), along with compounds associated with the projectile and gun barrel. The residue particles typically range in size from 1  $\mu\text{m}$  to over 100  $\mu\text{m}$  across, and can be examined easily with the Renishaw laboratory analyzer.



### Identification



Spectra of three gunshot residue particles, exhibiting features associated with barium oxide (BaO) and barium carbonate ( $\text{BaCO}_3$ ).

*Data courtesy of Prof. J. Yarwood, Sheffield Hallam University.*

# Forensic analysis

## Spectral libraries

The complexity of a molecule's Raman spectrum is directly related to the complexity of the molecule. Many molecules of interest to forensic scientists contain large numbers of atoms, and therefore have relatively complex spectra.

This richness allows materials to be distinguished by their spectra, but sometimes makes identification difficult. Spectral libraries are therefore a valuable tool, especially as they allow technicians (rather than Raman experts) to perform routine analyses.


### Available libraries

Renishaw has produced a variety of spectral libraries, each containing a large number of spectra. Those of particular relevance to forensic applications are:

- explosives
- narcotics
- inorganics (minerals, gems, oxides)
- polymers



### Identification of gems



Renishaw's spectral libraries are easy to use:

- acquire the spectrum and start the library search
- the results are ranked in order of quality of match, with the best first
- both the unidentified spectrum and the match are displayed on the computer
- the user can display other possible matches to confirm the computer's selection
- the match can be subtracted from the unknown to reveal other, weaker, constituent spectra
- the spectra and match statistics can be printed for later analysis

Spectrum from a supposed diamond is identified in seconds by the database as being that from a zircon.

The Renishaw inorganics database contains spectra from over 800 minerals, gemstones, and oxides.

## The Renishaw approach

Renishaw has long been world-renowned for its high precision measurement systems for quality control, research, and in-line process monitoring.

Renishaw was the first company to offer the new generation of Raman spectrometers, incorporating novel optical designs and recent advances in detector and filter technology. Renishaw's Raman systems are highly efficient, high throughput, compact, and easy to use instruments. They are the first choice for the Raman spectroscopist and analytical chemist alike.

These systems are causing a resurgence of interest in Raman spectroscopy, and are being applied routinely to both traditional and new application areas.

Renishaw continues to develop instruments for new applications, such as the ultraviolet Raman system for the analysis of **soils**, **contaminants**, and other difficult samples; please contact us for further details of our new systems.

Renishaw's Raman systems are currently in use with several law enforcement agencies around the world.

Renishaw is an ISO 9001 certified manufacturing company. Renishaw has been presented with a wide range of awards for technological and export achievement, both national and international.

## Global support & service

Renishaw recognises the importance of training, after sales service, and applications support, and is very proud of its world-wide servicing and product support network. Our team of Raman specialists and engineers offer rapid and comprehensive support:

- sample analysis
- "hands-on" technical training
- method development
- seminars and lectures, both on- and off-site



## Patents

Renishaw uses a variety of proprietary techniques to achieve the high instrument performance described in this brochure. These techniques are covered by a range of patents and patent applications:

|                 |                  |
|-----------------|------------------|
| CN 1162993      | US 5,194,912     |
| EP 0404901      | US 5,442,438     |
| EP 0542962      | US 5,689,333     |
| EP 0543578      | US 5,510,894     |
| EP 0638788      | US 5,818,047     |
| EP 0740133      | RU 97104125/4183 |
| EP 0776472      | WO 92/22793      |
| JP 2,733,491    | ZA 92/04145      |
| JP 2,834,954    | ZA 92/08761      |
| JP 504,650/1998 |                  |
| JP 119,865/1997 |                  |



# Global solutions from Renishaw

Renishaw is an established world leader in metrology, providing high performance, cost-effective solutions for measurement and increased productivity. A worldwide network of subsidiary companies and distributors provides exceptional service and support for its customers.

Renishaw designs, develops and manufactures products which conform to ISO 9001 standards.

Renishaw provides innovative solutions with the following products:

- Probe systems for inspection on CMMs (co-ordinate measuring machines).
- Systems for job set-up, tool setting and inspection on machine tools.
- Scanning and digitising systems.
- Laser and automated ballbar systems for performance measurement and calibration of machines.
- Encoder systems for high accuracy position feedback.
- Spectroscopy systems for non-destructive material analysis in laboratory and process environments.



Authorised distributor



**Renishaw plc**, Spectroscopy Products Division,  
Old Town, Wotton-under-Edge, Gloucestershire, GL12 7DW,  
United Kingdom  
Tel: +44 (0) 1453 844302 Fax +44 (0) 1453 844236



**Renishaw Inc.** 5277 Trillium Blvd., Hoffman Estates,  
IL 60192, USA  
Tel +1 847 843 3666  
Fax +1 847 843 1744



**Renishaw K. K.** Across City Nakano-Sakaue,  
38-1, Chuo 1-chome, Nakano-ku,  
Tokyo 164-0011, Japan  
Tel +81 3 5332 6021  
Fax +81 3 5332 6025



**Renishaw GmbH** Karl-Benz Straße 12,  
D 72124 Pliezhausen, Germany  
Tel +49 7127 9810  
Fax +49 7127 88237



**Renishaw S.A.** 15 rue Albert Einstein,  
Champs sur Marne,  
77437 Marne la Vallée, Cedex 2, France  
Tel +33 1 64 61 84 84  
Fax +33 1 64 61 65 26



**Renishaw S.p.A.** Via dei Prati 5,  
10044 Pianezza,  
Torino, Italy  
Tel +39 011 966 10 52  
Fax +39 011 966 40 83



**Renishaw Iberica S.A.** Edificio Océano,  
Calle Garrotxa 10-12,  
Parque Más Blau, 08820 Prat de Llobregat,  
Barcelona, Spain  
Tel +34 93 478 21 31  
Fax +34 93 478 16 08



**Renishaw A.G.** Poststrasse 5,  
CH 8808 Pfäffikon, Switzerland  
Tel +41 55 410 66 66  
Fax +41 55 410 66 69



**Renishaw (Hong Kong) Ltd.** Unit 4A,  
3/F, New Bright Building,  
11 Sheung Yuet Road, Kowloon Bay, Hong Kong  
Tel +852 2753 0638  
Fax +852 2756 8786



**Renishaw Latino Americana Ltda.**  
Calçada dos Crisântemos 22,  
Centro Comercial de Alphaville,  
C.e.p. 06453-000, Barueri-SP, Brazil  
Tel +55 11 7295 2866  
Fax +55 11 7295 1641

## Renishaw's Representative Offices



**The Peoples' Republic of China**, Beijing  
Tel +86 10 8448 5306  
Fax +86 10 84481528



**The Peoples' Republic of China**, Shanghai  
Tel +86 21 6218 5183  
Fax +86 21 6218 5193



**Singapore**  
Tel +65 8975466  
Fax +65 8975467



**Indonesia**, Jakarta  
Tel +62 21 428 70153  
Fax +62 21 424 3934



**Australia**, Melbourne  
Tel +61 3 9553 8267  
Fax +61 3 9592 6738



## Renishaw's Liaison Offices



**India**, Bangalore  
Tel +91 80 5095419  
Fax +91 80 5095421



**South Korea**, Seoul  
Tel +82 2 565 6878  
Fax +82 2 565 6879

# RENISHAW

email: [raman@renishaw.com](mailto:raman@renishaw.com)

Internet: <http://www.renishaw.com>